

Ko Validated Anti-Jak2 Rabbit Polyclonal Antibody For Wb, Elisa - O60674

Item Number: CM0030026



Introduction

KO validated rabbit polyclonal antibody targeting human JAK2, suitable for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for signal transduction, chromatin regulation, and cytokine signaling studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] JAK2 Rabbit pAb
Remarks/Alias	JTK10; JAK2
Species	Human
Gene ID (Human)	3717
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 933-1132 of human JAK2 (NP_004963.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60674
Protein Weight	131 kDa
Shipping	Ice bag