

[KO Validated] Anti-Keap1 Rabbit Monoclonal Antibody (Inrf2/Klhl19) - Q14145

Item Number: CM0009101



Introduction

Validated KEAP1 (INrf2/KLHL19) rabbit monoclonal antibody for WB, IF/ICC, ELISA. Suitable for human, mouse, rat. Key oxidative stress sensor; E3 ubiquitin ligase adapter for NRF2 degradation.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] KEAP1 Rabbit mAb
Remarks/Alias	INrf2, KLHL19
Species	Human
GeneID (Human)	9817
GeneID	9817
Immunogen	Synthetic peptide within amino acids 525-624 of human KEAP1 (NP_036421.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14145
Protein Weight	70kDa
Shipping	Ice bag