

[KO Validated] Anti-Ki67 Monoclonal Antibody For Wb, If/Icc, Ihc-P, Mihc, Elisa - P46013

Item Number: CM0002758



Introduction

Rabbit monoclonal antibody targeting human Ki67, a proliferation marker. Validated for WB, IF/ICC, IHC-P, mIHC, and ELISA. Cross-reacts with mouse and rat. Suitable for cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Ki67 Rabbit mAb
Remarks/Alias	KiA; MIB-; MIB-1; PPP1R105; Ki67
Species	Human
GeneID (Human)	4288
GeneID	4288
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1200-1300 of human Ki67 (NP_002408.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, mIHC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46013
Protein Weight	359 kDa
Shipping	Ice bag