

Ko Validated Anti-Kmt5C Polyclonal Antibody For Wb, Elisa - Q86Y97

Item Number: CM0023237



Introduction

KO validated anti-KMT5C polyclonal antibody validated for WB and ELISA. Reacts with human and mouse. Targets histone methyltransferase KMT5C (SUV420H2) for chromatin regulation and DNA repair studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] KMT5C Rabbit pAb
Remarks/Alias	SUV420H2; Suv4-20h2; 5C
Species	Human
GeneID (human)	84787
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 100-200 of human KMT5C (NP_116090.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q86Y97
Protein Weight	52kDa
Shipping	Ice bag