

[KO Validated] Anti-Kras Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P01116

Item Number: CM0006359



Introduction

KO validated anti-KRAS rabbit monoclonal antibody for WB, IHC-P, and ELISA. Detects human, mouse, and rat KRAS, a key Ras-MAPK signaling protein involved in cell proliferation and oncogenesis.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] KRAS Rabbit mAb
Remarks/Alias	NS; NS3; OES; CFC2; RALD; K-Ras; KRAS1; KRAS2; RASK2; KI-RAS; C-K-RAS; K-RAS2A; K-RAS2B; K-RAS4A; K-RAS4B; K-Ras 2; 'C-K-RAS; c-Ki-ras; c-Ki-ras2; AS
Species	Human
GeneID (Human)	3845
GeneID	3845
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01116
Protein Weight	22kDa
Shipping	Ice bag