

KO Validated Anti-Kras Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P01116

Item Number: CM0017525



Introduction

KO Validated KRAS Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets GTPase KRas, a key signal transducer in the Ras-MAPK pathway. Suitable for cancer and cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] KRAS Rabbit pAb
Remarks/Alias	NS; NS3; OES; CFC2; RALD; K-Ras; KRAS1; KRAS2; RASK2; KI-RAS; C-K-RAS; K-RAS2A; K-RAS2B; K-RAS4A; K-RAS4B; K-Ras 2; 'C-K-RAS; c-Ki-ras; c-Ki-ras2; AS
Species	Human
GeneID (human)	3845
GeneID	3845
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-188 of human KRAS (NP_203524.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01116
Protein Weight	22kDa
Shipping	Ice bag