


Ko Validated Anti-Lactate Dehydrogenase B/Ldh-B Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P07195

Item Number: CM0008552



Introduction

KO Validated Lactate Dehydrogenase B/LDH-B rabbit monoclonal antibody for WB, IF/ICC, IF-P, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for cardiac muscle glycolysis studies. CamelBio supplied.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Lactate Dehydrogenase B/LDH-B Rabbit mAb
Remarks/Alias	LDH-B; LDH-H; LDHBD; TRG-5; HEL-S-281
Species	Human
Gene ID	3945
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2-93 of human LDHB (NP_002291.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07195
Protein Weight	37kDa
Shipping	Ice bag