

Ko-Validated Anti-Lamin A/C Mouse Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P02545

Item Number: CM0008599



Introduction

KO-validated mouse monoclonal anti-Lamin A/C antibody for WB, IHC-P, IF/ICC, ELISA. Detects lamin A/C in human, mouse, and rat samples. Useful for nuclear lamina, laminopathy, and cellular senescence studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Lamin A/C Mouse mAb
Remarks/Alias	FPL; IDC; LFP; CDDC; EMD2; FPLD; HGPS; LDP1; LMN1; LMNC; MADA; PRO1; CDCD1; CMD1A; FPLD2; LMNL1; CMT2B1; LGMD1B
Species	Human
GeneID	4000
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 403-572 of human Lamin A/C (NP_733821.1).
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02545
Protein Weight	74kDa
Shipping	Ice bag