

[KO Validated] Anti-Lamin A/C Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P02545

Item Number: CM0024355



Introduction

KO-validated rabbit polyclonal antibody against Lamin A/C (P02545). Validated for WB, IHC-P, IF/ICC, IP, ELISA. Reacts with human, mouse, rat. Ideal for nuclear lamina and laminopathy research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Lamin A/C Rabbit pAb
Remarks/Alias	FPL; IDC; LFP; CDDC; EMD2; FPLD; HGPS; LDP1; LMN1; LMNC; MADA; PRO1; CDCD1; CMD1A; FPLD2; LMNL1; CMT2B1; LGMD1B; /C
Species	Human
GeneID (Human)	4000
GeneID	4000
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 403-572 of human Lamin A/C (NP_733821.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02545
Protein Weight	74kDa
Shipping	Ice bag