

## [KO Validated] Anti-Lamp2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P13473

Item Number: CM0003230



### Introduction

KO validated anti-LAMP2 rabbit monoclonal antibody for WB, IHC-P, ELISA. Detects human, mouse, rat LAMP2. Suitable for autophagy, lysosome biogenesis, and Danon disease research. UniProt P13473.

[Learn More](#)

| Parameter        | Value                                                                                               |
|------------------|-----------------------------------------------------------------------------------------------------|
| Product Name     | [KO Validated] LAMP2 Rabbit mAb                                                                     |
| Remarks/Alias    | DND; LAMPB; CD107b; LAMP-2; LGP-96; LGP110; LAMP2                                                   |
| Species          | Human                                                                                               |
| GeneID (human)   | 3920                                                                                                |
| GeneID           | 3920                                                                                                |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 311-410 of human LAMP2 (P13473). |
| Source           | Rabbit                                                                                              |
| Category         | Monoclonal Antibodies                                                                               |
| Application      | WB, IHC-P, ELISA                                                                                    |
| Cross Reactivity | Human, Mouse, Rat                                                                                   |
| SWISS            | P13473                                                                                              |
| Protein Weight   | 45kDa                                                                                               |
| Shipping         | Ice bag                                                                                             |