

## [KO Validated] Anti-Lamp2 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Fc, Elisa - P13473

Item Number: CM0012170



### Introduction

High-specificity Anti-LAMP2 Rabbit Monoclonal Antibody, KO validated, suitable for WB, IHC-P, IF/ICC, FC and ELISA. Ideal for studying lysosome biogenesis, autophagy, and lysosomal storage diseases.

[Learn More](#)

| Parameter        | Value                                                                                                                                    |
|------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | [KO Validated] LAMP2 Rabbit mAb                                                                                                          |
| Remarks/Alias    | DND; LAMPB; CD107b; LAMP-2; LGP-96; LGP110                                                                                               |
| Species          | Human                                                                                                                                    |
| GeneID (Human)   | 3920                                                                                                                                     |
| GeneID           | 3920                                                                                                                                     |
| Immunogen        | Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 29-375 of human LAMP2 (NP_002285.1). |
| Source           | Rabbit                                                                                                                                   |
| Category         | Monoclonal Antibodies                                                                                                                    |
| Application      | WB, IHC-P, IF/ICC, FC (intra), ELISA                                                                                                     |
| Cross Reactivity | Human                                                                                                                                    |
| SWISS            | P13473                                                                                                                                   |
| Protein Weight   | 45kDa                                                                                                                                    |
| Shipping         | Ice bag                                                                                                                                  |