


[KO Validated] Anti-Ldhb Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P07195

Item Number: CM0022351



Introduction

KO-validated rabbit polyclonal antibody targeting human LDHB (Lactate Dehydrogenase B). Validated for WB, IF-P, ELISA in human, mouse, rat. LDHB catalyzes pyruvate-lactate interconversion and is highly expressed in cardiac muscle.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] LDHB Rabbit pAb
Remarks/Alias	LDH-B; LDH-H; LDHBD; TRG-5; HEL-S-281; HB
Species	Human
Gene ID (Human)	3945
Gene ID	3945
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 2-334 of human LDHB (NP_002291.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07195
Protein Weight	37kDa
Shipping	Ice bag