

## [Ko Validated] Anti-Litaf Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q99732

Item Number: CM0027160

### Introduction

KO-validated rabbit polyclonal antibody against human LITAF, suitable for WB, IHC-P, IP, and ELISA. Detects 17 kDa protein in human and mouse samples. Ideal for endosomal trafficking and cytokine research.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | [KO Validated] LITAF Rabbit pAb                                                                                  |
| Remarks/Alias    | PIG7; SIMPLE; TP53I7; AF                                                                                         |
| Species          | Human                                                                                                            |
| Gene ID (Human)  | 9516                                                                                                             |
| Gene ID          | 9516                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-161 of human LITAF (NP_004853.2) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, IHC-P, IP, ELISA                                                                                             |
| Cross Reactivity | Human, Mouse                                                                                                     |
| SWISS            | Q99732                                                                                                           |
| Protein Weight   | 17kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |