


[KO Validated] Anti-Lsh/Hells Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q9Nrz9

Item Number: CM0025089



Introduction

KO-validated rabbit polyclonal antibody targeting human Lsh/HELLS (Q9NRZ9). Suitable for WB, IP, and ELISA with cross-reactivity to mouse and rat. An ideal tool for chromatin remodeling, DNA methylation, and genomic stability studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Lsh/HELLS Rabbit pAb
Remarks/Alias	LSH; ICF4; PASG; SMARCA6; Nbla10143; LS
Species	Human
Gene ID (Human)	3070
Gene ID	3070
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 719-838 of human Lsh/HELLS (NP_060533.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NRZ9
Protein Weight	97 kDa
Shipping	Ice bag