


Ko Validated Anti-MacroH2A.1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O75367

Item Number: CM0029408



Introduction

Rabbit polyclonal antibody against human macroH2A.1 (KO validated), suitable for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. For research use in chromatin biology and epigenetics.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] macroH2A.1 Rabbit pAb
Remarks/Alias	H2A.y; H2A/y; H2AFY; mH2A1; H2AF12M; MACROH2A1.1; macroH2A1.2; .1
Species	Human
GeneID (Human)	9555
GeneID	9555
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 123-372 of human macroH2A.1 (NP_613258.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75367
Protein Weight	39kDa
Shipping	Ice bag