

[KO Validated] Anti-Mapkapk3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q16644

Item Number: CM0021577



Introduction

Rabbit polyclonal antibody against MAPKAPK3, validated for WB, IF/ICC, ELISA; reacts with human, mouse, rat. Stress-activated kinase involved in cytokine production, endocytosis, and chromatin remodeling.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] MAPKAPK3 Rabbit pAb
Remarks/Alias	3PK; MK3; MK-3; MDPT3; MAPKAP3; MAPKAP-K3; MAPKAPK-3; K3
Species	Human
GeneID (human)	7867
GeneID	7867
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human MAPKAPK3 (NP_004626.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16644
Protein Weight	43kDa
Shipping	Ice bag