


Ko Validated Anti-MLh1 Monoclonal Antibody For Wb, Ihc-P, Elisa - P40692

Item Number: CM0007468



Introduction

KO validated rabbit monoclonal antibody targeting human MLH1, suitable for WB, IHC-P, and ELISA applications. Cross-reacts with mouse. Ideal for DNA mismatch repair and cancer research studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] MLH1 Rabbit mAb
Remarks/Alias	FCC2; COCA2; HNPCC; MLH-1; hMLH1; HNPCC2; LYNCH2; MMRCS1
Species	Human
GeneID (Human)	4292
GeneID	4292
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 350-450 of human MLH1 (NP_000240.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P40692
Protein Weight	85 kDa
Shipping	Ice bag