


[Ko Validated] Anti-Msh2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P43246

Item Number: CM0005473



Introduction

KO validated anti-MSH2 rabbit monoclonal antibody for WB, IHC-P, ELISA. Reacts with human, mouse, rat MSH2. DNA mismatch repair protein Msh2, for cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] MSH2 Rabbit mAb
Remarks/Alias	FCC1; COCA1; HNPCC; LCF52; MSH-2; hMSH2; HNPCC1; LYNCH1; MMRCS2; H2
Species	Human
GeneID	4436
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 330-583 of human MSH2 (NP_000242.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43246
Protein Weight	105kDa
Shipping	Ice bag