

[Ko Validated] Anti-Mta2 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa, Chip - O94776

Item Number: CM0008627

Introduction

[KO Validated] Rabbit monoclonal antibody targeting human MTA2, a component of the NuRD complex involved in chromatin remodeling. Validated for WB, IHC-P, IF/ICC, ELISA, and CHIP.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] MTA2 Rabbit mAb
Remarks/Alias	PID; MTA1L1; A2
Species	Human
GeneID (Human)	9219
GeneID	9219
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 569-668 of human MTA2 (O94776)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA, CHIP
Cross Reactivity	Human
SWISS	O94776
Protein Weight	75kDa
Shipping	Ice bag