

Ko Validated Anti-Mta2 Rabbit Polyclonal Antibody For Wb, Elisa, Chip - O94776

Item Number: CM0019217

Introduction

KO Validated Anti-MTA2 rabbit polyclonal antibody, validated for Western blot, ELISA, and ChIP. Detects MTA2 in human, mouse, and rat samples. Suitable for chromatin remodeling and transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] MTA2 Rabbit pAb
Remarks/Alias	PID; MTA1L1; A2
Species	Human
GeneID (human)	9219
GeneID	9219
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-180 of human MTA2 (NP_004730.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	O94776
Protein Weight	75kDa
Shipping	Ice bag