


[KO Validated] Anti-Mth1 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P36639

Item Number: CM0016677



Introduction

Rabbit polyclonal antibody targeting human MTH1/NUDT1 (oxidized purine nucleoside triphosphate hydrolase). Validated for knockout studies, suitable for WB, IP, ELISA applications. Detects 18kDa protein, cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] MTH1 Rabbit pAb
Remarks/Alias	MTH1; H1
Species	Human
GeneID (Human)	4521
GeneID	4521
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-179 of human MTH1 (NP_945192.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	P36639
Protein Weight	18kDa
Shipping	Ice bag