

[Ko Validated] Anti-N-Cadherin Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P19022

Item Number: CM0002400

Introduction

Rabbit monoclonal antibody validated for knockout applications against N-Cadherin (CDH2). Suitable for WB, IF-P, IHC-P, and ELISA with cross-reactivity to human, mouse, and rat. Ideal for cell adhesion and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] N-Cadherin Rabbit mAb
Remarks/Alias	CDHN, NCAD, ACOGS, ADHD8, CD325, ARVD14, CDw325
Species	Human
Gene ID	1000
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 100-200 of human N Cadherin (P19022).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19022
Protein Weight	100 kDa
Shipping	Ice bag