

Ko Validated Anti-N-Terminal Alpha-Synuclein Monoclonal Antibody [Mjff-B-3A4] For Wb, If/Icc, Elisa - P37840

Item Number: CM0010732



Introduction

KO validated rabbit monoclonal antibody targeting N-terminal alpha-synuclein (SNCA/PARK1). Validated for WB, IF/ICC, ELISA; cross-reacts with human, mouse, rat. Ideal for synaptic function and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] N-terminal alpha-synuclein mAb [MJFF-B-3A4]
Remarks/Alias	NACP; PARK1; α -Synuclein; SNCA
Species	Human
Gene ID	6622
Immunogen	A synthetic peptide corresponding to amino acids 40-50 of human alpha-synuclein (NP_000336.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P37840
Protein Weight	14kDa
Shipping	Ice bag