


[KO Validated] Anti-Ndufc2 Polyclonal Antibody For Wb And Elisa - O95298

Item Number: CM0018032



Introduction

High-quality rabbit polyclonal antibody against human NDUFC2, a key complex I subunit. Validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for mitochondrial research and respiratory chain studies.

[Learn More](#)

Parameter	Specification
Product Name	[KO Validated] NDUFC2 Rabbit pAb
Remarks/Alias	HLC-1; B14.5b; MC1DN36; NADHDH2; CI-B14.5b; C2
Species	Human
Gene ID (Human)	4718
Gene ID	4718
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human NDUFC2 (NP_004540.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95298
Protein Weight	14kDa
Shipping	Ice bag