


[KO Validated] Anti-Nosip Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y314

Item Number: CM0022895



Introduction

Rabbit polyclonal antibody against human NOSIP (nitric oxide synthase-interacting protein, E3 ubiquitin ligase), validated for WB and ELISA applications. Suitable for research on nitric oxide regulation and development. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] NOSIP Rabbit pAb
Remarks/Alias	CGI-25; NOSIP
Species	Human
GeneID (Human)	51070
GeneID	51070
Immunogen	Recombinant protein corresponding to a sequence within amino acids 1-301 of human NOSIP (NP_057037.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y314
Protein Weight	33kDa
Shipping	Ice bag