

[Ko Validated] Anti-Nrd1 Rabbit Polyclonal Antibody For Wb And Elisa - O43847

Item Number: CM0021379

Introduction

KO-validated anti-NRD1 rabbit polyclonal antibody for WB and ELISA. Detects nardilysin (~132 kDa) in human, mouse, and rat samples. Ideal for axonal maturation and myelination studies. UniProt O43847.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] NRD1 Rabbit pAb
Remarks/Alias	NRD1; hNRD1; hNRD2; D1
Species	Human
GeneID (Human)	4898
GeneID	4898
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 160-400 of human NRD1 (NP_001095132.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43847
Protein Weight	132kDa
Shipping	Ice bag