


[KO Validated] Anti-Nsdhl Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q15738

Item Number: CM0020012



Introduction

KO validated rabbit polyclonal antibody targeting human NSDHL (sterol-4-alpha-carboxylate 3-dehydrogenase), involved in cholesterol biosynthesis. Suitable for WB, IHC-P, IF/ICC, ELISA; reacts with human, mouse, rat. - Q15738

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] NSDHL Rabbit pAb
Remarks/Alias	H105E3; XAP104; SDR31E1; HL
Species	Human
GeneID (Human)	50814
GeneID	50814
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human NSDHL (NP_057006.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15738
Protein Weight	42kDa
Shipping	Ice bag