

[KO Validated] Anti-Nucleobindin-1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q02818

Item Number: CM0026414



Introduction

KO-validated rabbit polyclonal antibody targeting nucleobindin-1 (NUCB1, CALNUC), suitable for WB, IF/ICC, ELISA. Cross-reactivity with human, mouse, rat. Validated for research on Golgi calcium homeostasis.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] NUCB1 Rabbit pAb
Remarks/Alias	NUC; CALNUC; B1
Species	Human
Gene ID (Human)	4924
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 302-461 of human NUCB1 (NP_006175.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02818
Protein Weight	54kDa
Shipping	Ice bag