

[KO Validated] Anti-Pak2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q13177

Item Number: CM0023891



Introduction

High-quality rabbit polyclonal antibody against human PAK2, validated for WB, IF/ICC, and ELISA. Recognizes mouse and rat PAK2. Ideal for studying cytoskeleton regulation, apoptosis, and cell signaling.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PAK2 Rabbit pAb
Remarks/Alias	KNO2; PAK65; PAKgamma; K2
Species	Human
GeneID (Human)	5062
GeneID	5062
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-212 of human PAK2 (NP_000903.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13177
Protein Weight	58kDa
Shipping	Ice bag