


Ko Validated Anti-Pdcd6 Polyclonal Antibody (Alg2) For Wb, Ip, Elisa - O75340

Item Number: CM0025571



Introduction

Rabbit polyclonal antibody against PD CD6 (ALG-2), validated for Western blot, immunoprecipitation, and ELISA. Cross-reacts with human, mouse, and rat. Suitable for apoptosis and calcium signaling studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PD CD6 Rabbit pAb
Remarks/Alias	ALG2; ALG-2; PEF1B; D6
Species	Human
Gene ID (Human)	10016
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-191 of human PD CD6 (NP_037364.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75340
Protein Weight	22kDa
Shipping	Ice bag