

[Ko Validated] Anti-Pepck/Pck2 Rabbit Polyclonal Antibody - Q16822

Item Number: CM0020733

Introduction

KO validated rabbit polyclonal antibody targeting human PEPCK/PCK2. Reacts with human, mouse, rat. Suitable for WB, IHC-P, IF/ICC, IP, ELISA. Research tool for gluconeogenesis, ferroptosis, and mitochondrial metabolism.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PEPCK/PCK2 Rabbit pAb
Remarks/Alias	PEPCK; PEPCK2; PEPCK-M; K2
Species	Human
GeneID (Human)	5106
GeneID	5106
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human PEPCK/PCK2 (NP_004554.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16822
Protein Weight	71kDa
Shipping	Ice bag