

Ko Validated Anti-Phd2/Egln1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Gzt9

Item Number: CM0018123



Introduction

Rabbit polyclonal antibody targeting human PHD2/EGLN1, KO validated. Suitable for WB, IHC-P, IF/ICC, and ELISA. Recognizes human, mouse, and rat PHD2. Ideal for hypoxia and HIF signaling research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PHD2/EGLN1 Rabbit pAb
Remarks/Alias	HPH2; PHD2; SM20; ECT3; HALAH; HPH-2; HIFPH2; ZMYND6; C1orf12; HIF-PH2; N1
Species	Human
GeneID (human)	54583
GeneID	54583
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZT9
Protein Weight	46kDa
Shipping	Ice bag