

[KO Validated] Anti-Pka Rii α (Prkar2A) Polyclonal Antibody - P13861

Item Number: CM0015737



Introduction

[KO Validated] Rabbit polyclonal antibody against human PKA RII α (PRKAR2A). Suitable for Western blot, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for studying cAMP signaling and PKA regulation.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PKA RII α (PRKAR2A) Rabbit pAb
Remarks/Alias	PKR2; PRKAR2; A)
Species	Human
Gene ID (Human)	5576
Gene ID	5576
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-404 of human PKA RII α (PRKAR2A)/PKR2 (NP_004148.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13861
Protein Weight	46kDa
Shipping	Ice bag