

[KO Validated] Anti-Pml Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - P29590

Item Number: CM0011870



Introduction

Rabbit monoclonal antibody against human PML protein, validated for knockout applications. Suitable for WB, IP, IF/ICC, IHC-P, and ELISA. Recognizes human PML, relevant in tumor suppression, apoptosis, and antiviral defense.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PML Rabbit mAb
Remarks/Alias	MYL; RNF71; PP8675; TRIM19
Species	Human
GeneID (human)	5371
GeneID	5371
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-185 of human PML (NP_150241.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P29590
Protein Weight	98kDa
Shipping	Ice bag