

Ko Validated Anti-Pnp Rabbit Polyclonal Antibody For Wb, Elisa - P00491

Item Number: CM0023608



Introduction

KO-validated rabbit polyclonal antibody against human PNP (purine nucleoside phosphorylase) for WB and ELISA. Cross-reacts with mouse. Molecular weight ~32kDa. Suitable for nucleoside metabolism and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PNP Rabbit pAb
Remarks/Alias	NP; PUNP; PRO1837
Species	Human
GeneID (Human)	4860
GeneID	4860
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-289 of human PNP (NP_000261.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P00491
Protein Weight	32kDa
Shipping	Ice bag