


KO Validated Anti-Prmt4/Carm1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q86X55

Item Number: CM0022784



Introduction

KO validated rabbit polyclonal antibody targeting PRMT4/CARM1 (SWISS: Q86X55), suitable for WB, IHC-P, and ELISA. Detects human, mouse, and rat proteins; ideal for transcription and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PRMT4/CARM1 Rabbit pAb
Remarks/Alias	PRMT4; M1
Species	Human
Gene ID (Human)	10498
Gene ID	10498
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 486-603 of human PRMT4/CARM1 (NP_954592.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86X55
Protein Weight	66kDa
Shipping	Ice bag