


[KO Validated] Anti-Psap Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P07602

Item Number: CM0021050



Introduction

Rabbit polyclonal antibody against PSAP (Prosaposin), KO validated, for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Essential for lysosomal storage disease and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PSAP Rabbit pAb
Remarks/Alias	GLBA; SAP1; SAP2; PSAPD; PARK24; AP
Species	Human
GeneID	5660
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB,IHC-P, IF/ICC, ELISA
Cross Reactivity	Human,Mouse,Rat
SWISS	P07602
Protein Weight	58kDa
Shipping	Ice bag