

## [KO Validated] Anti-Ptbp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P26599

Item Number: CM0024224



### Introduction

KO validated PTBP1 rabbit polyclonal antibody for western blot, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat PTBP1, a key splicing regulator. Ideal for IVD and research applications.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | [KO Validated] PTBP1 Rabbit pAb                                                                                  |
| Remarks/Alias    | PTB; PTB2; PTB3; PTB4; pPTB; HNRPI; PTB-1; PTB-T; HNRNPI; HNRNP-I; P1                                            |
| Species          | Human                                                                                                            |
| GeneID (Human)   | 5725                                                                                                             |
| GeneID           | 5725                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human PTBP1 (NP_114368.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | P26599                                                                                                           |
| Protein Weight   | 60kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |