

[KO Validated] Anti-Rad23A Rabbit Polyclonal Antibody For Wb, Elisa - P54725

Item Number: CM0020598



Introduction

KO validated rabbit polyclonal antibody against human Rad23A. Suitable for WB and ELISA applications, with cross-reactivity to human and mouse. Ideal for DNA repair and demethylation research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Rad23A Rabbit pAb
Remarks/Alias	HR23A; HHR23A; 3A
Species	Human
GeneID (Human)	5886
GeneID	5886
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-363 of human Rad23A (NP_005044.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P54725
Protein Weight	40kDa
Shipping	Ice bag