

[KO Validated] Anti-Rbm3 Polyclonal Antibody For Wb, If/Icc, Elisa - P98179

Item Number: CM0021767



Introduction

RBM3 rabbit polyclonal antibody validated for Western blot, IF/ICC, and ELISA. Reacts with human, mouse, rat. Targets cold-inducible RNA-binding protein 3 (P98179) involved in translational regulation and stress response. Suitable for research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] RBM3 Rabbit pAb
Remarks/Alias	RNPL; IS1-RNPL; M3
Species	Human
GeneID (Human)	5935
GeneID	5935
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-157 of human RBM3 (NP_006734.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P98179
Protein Weight	17kDa
Shipping	Ice bag