

[KO Validated] Anti-Rela Rabbit Polyclonal Antibody For Wb, Elisa - Q04207

Item Number: CM0022670



Introduction

Rabbit polyclonal antibody against mouse RELA (NF-kB p65), validated for Western blot and ELISA, cross-reacts with human. Ideal for studying inflammation, immunity, cell growth, and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] RELA Rabbit pAb
Remarks/Alias	p65, p65 NFkB, p65 NF-kappa B, LA
Species	Mouse
GeneID	19697 (Mouse)
Immunogen	Recombinant protein of mouse RELA
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q04207
Protein Weight	60kDa
Shipping	Ice bag