

[Ko Validated] Anti-Rhogdi Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P52565

Item Number: CM0015656



Introduction

Polyclonal rabbit antibody against human RhoGDI (ARHGDIA). Validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Study Rho protein regulation and cell motility in cancer and normal cells.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] RhoGDI Rabbit pAb
Remarks/Alias	GDIA1; NPHS8; RHOGDI; RHOGDI-1; HEL-S-47e; DI
Species	Human
GeneID (Human)	396
GeneID	396
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-204 of human RhoGDI (NP_001172006.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52565
Protein Weight	23kDa
Shipping	Ice bag