


KO Validated Anti-Rpl22 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P35268

Item Number: CM0026013



Introduction

KO validated rabbit polyclonal antibody against ribosomal protein L22 (RPL22), for Western blot, IHC-P, and ELISA. Cross-reacts with human, mouse, rat. Suitable for ribosome biology and translational control research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] RPL22 Rabbit pAb
Remarks/Alias	EAP; L22; HBP15; HBP15/L22; 22
Species	Human
GeneID (Human)	6146
GeneID	6146
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-122 of human RPL22 (NP_000974.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35268
Protein Weight	15kDa
Shipping	Ice bag