

[Ko Validated] Anti-Scyl1 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Kg9

Item Number: CM0026083

Introduction

High-quality rabbit polyclonal antibody targeting SCYL1 (N-terminal kinase-like protein). Validated for Western blot and ELISA across human, mouse, and rat samples. Suitable for trafficking and transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SCYL1 Rabbit pAb
Remarks/Alias	GKLP; NKTl; NTKL; P105; TAPK; TEIF; TRAP; HT019; SCAR21; L1
Species	Human
GeneID (Human)	57410
GeneID	57410
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 460-740 of human SCYL1 (NP_065731.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96KG9
Protein Weight	90kDa
Shipping	Ice bag