


[KO Validated] Anti-Sec61B Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P60468

Item Number: CM0013841



Introduction

Validated rabbit polyclonal antibody targeting human SEC61B, a core subunit of the ER translocon complex. Suitable for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for ER translocation and protein trafficking studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SEC61B Rabbit pAb
Remarks/Alias	SEC61B; 1B
Species	Human
GeneID	10952
Immunogen	Synthetic peptide within amino acids 1-96 of human SEC61B (NP_006799.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60468
Protein Weight	10kDa
Shipping	Ice bag