

[Ko Validated] Anti-Selenbp1 Rabbit Polyclonal Antibody - Q13228

Item Number: CM0013506

Introduction

[KO Validated] Rabbit polyclonal antibody against SELENBP1 (methanethiol oxidase) for WB, IF/ICC, ELISA. Validated in human and mouse samples. Ideal for cancer and neuroscience research.

[Learn More](#)

Parameter	Specification
Product Name	[KO Validated] SELENBP1 Rabbit pAb
Remarks/Alias	MTO; LPSB; SP56; hSBP; EHMT0; SBP56; HEL-S-134P; P1
Species	Human
GeneID (Human)	8991
GeneID	8991
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 243-472 of human SELENBP1 (NP_003935.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13228
Protein Weight	52kDa
Shipping	Ice bag