


[KO Validated] Anti-Sephs1 Rabbit Polyclonal Antibody - P49903

Item Number: CM0016713



Introduction

KO-validated anti-SEPHS1 rabbit polyclonal antibody for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat SEPHS1. Useful for zincore complex and transcription factor stability research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SEPHS1 Rabbit pAb
Remarks/Alias	SPS; SELD; SPS1; S1
Species	Human
GeneID (Human)	22929
GeneID	22929
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-392 of human SEPHS1 (NP_036379.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49903
Protein Weight	43kDa
Shipping	Ice bag