

[KO Validated] Anti-Serpinb5 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P36952

Item Number: CM0015315



Introduction

Rabbit polyclonal antibody targeting human SERPINB5 (maspin) validated by WB, IHC-P, ELISA. Cross-reacts with mouse. Supports tumor suppressor and cancer research studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SERPINB5 Rabbit pAb
Remarks/Alias	PI5; maspin; B5
Species	Human
Gene ID	5268
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-189 of human SERPINB5 (NP_002630.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P36952
Protein Weight	42kDa
Shipping	Ice bag