


[KO Validated] Anti-Setd7 Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Wts6

Item Number: CM0027458



Introduction

KO validated rabbit polyclonal antibody against SETD7, a histone H3K4 methyltransferase. Plays a key role in epigenetic transcriptional activation and methylation of non-histone substrates. Suitable for WB, IF/ICC, ELISA in human, mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SETD7 Rabbit pAb
Remarks/Alias	KMT7; SET7; SET9; SET7/9; D7
Species	Human
GeneID	80854
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-366 of human SETD7 (NP_085151.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WTS6
Protein Weight	41kDa
Shipping	Ice bag