

[Ko Validated] Anti-Sirt1 Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q96Eb6

Item Number: CM0005989

Introduction

KO-validated rabbit monoclonal antibody against human SIRT1 (SIR2/SIR2L1). Suitable for WB, IF/ICC, IHC-P, FC, and ELISA. Detects human, mouse, and rat SIRT1, with an ~82 kDa band. Ideal for epigenetics and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SIRT1 Rabbit mAb
Remarks/Alias	SIR2; SIR2L1; SIR2alpha
Species	Human
GeneID (Human)	23411
GeneID	23411
Immunogen	Recombinant protein of human SIRT1.
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, FC (intra), ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96EB6
Protein Weight	82 kDa
Shipping	Ice bag